



LOST YOUR REMOTE?

- 1** Remove the Paraglide device from the wheelchair. Open the green charge port cap and find the small button near the charging port. Press down on this button briefly (about 1 second) and the button will illuminate and flash. This means the device has sent a signal to the remote control.
- 2** Once the remote receives the signal it will sound an alarm that will continue for two minutes or until the "Reset" button is pressed. The remote can take up to 30 seconds to receive the signal.

DEVICE NOT FOUND?

If your remote displays the message "Device Not Found", first check to ensure your Paraglide device is charged by pressing the button beside the charge port. If it illuminates, showing your device is charged, follow these steps:

- 1** Press and hold the small button near the charging port for 3-5 seconds or until the blue LED light begins to flash.
- 2** Using the remote, go to the menu page by holding the Data and Reset buttons. Select the option Reconnect by pressing the Data button. Your remote will now reconnect to the device.
- 3** If unsuccessful, you'll need to perform a system reset. See pg. 39 in operators manual for instructions.